

Bioseal®

Sterile Cast Padding Products



Bioseal produces a full line of Sterile Cast Padding to fill the needs of the Orthopaedic Surgeon. This includes BSN's Sof-Rol and Specialist, Kendall's Webril line, and our own Bi-Web Cotton Cast Padding & Bioline Synthetic Cast Padding. All of our Cast Paddings are radiation sterilized to insure maximum patient safety and no risk of ETO residual.

Bi-Web

- Best value
- 100% cotton fiber for comfort
- Sterile or non-sterile

Webril™

- Sized finish for added strength
- 100% cotton fiber for comfort
- Kendall "Brand Name" quality

Bioline Synthetic

- Designed for use under new hi-tech casts

Sof-Rol™

- 100% rayon fiber
- "Brand Name" quality

Specialist™

- Pleated design for added conformability
- 100% cotton fiber for comfort
- "Brand Name" quality

*Webril is a trademark of The Kendall Company
Specialist and Sof-Rol are trademarks of BSN*



Sterile Cast Padding Products

Bi-Web

Sterile

3024/24	2" x 4 yds
3034/24	3" x 4 yds
3044/24	4" x 4 yds
3064/24	6" x 4 yds

Webril™

Sterile

3024k/24	2" x 4 yds
3034k/24	3" x 4 yds
3044k/24	4" x 4 yds
3064k/24	6" x 4 yds

Synthetic

Sterile

3534/20	3" x 4 yds	20/cs
3544/20	4" x 4 yds	20/cs
3564/20	6" x 4 yds	20/cs

Sof-Rol™

Sterile

24/cs	3124/24	2" x 4 yds	24/cs
24/cs	3134/24	3" x 4 yds	24/cs
24/cs	3144/24	4" x 4 yds	24/cs
24/cs	3164/24	6" x 4 yds	24/cs

Specialist™

Sterile

24/cs	3224/24	2" x 4 yds	24/cs
24/cs	3234/24	3" x 4 yds	24/cs
24/cs	3244/24	4" x 4 yds	24/cs
24/cs	3264/24	6" x 4 yds	24/cs